

MEDECO B.V.

Product Specifications

Brand : Klinion® Protection → POLYMERCOATED
Product : Latex, powder-free, Polymer Coated, non-sterile, examination gloves

Packaging & Product Reference

Product Reference	Description	Size	Color	Box Dimension (mm)*	Piece(s) per box	Outer Carton Dimension (mm)*	Box(es) per Outer Carton
102645	Klinion Personal Protection Latex Polymer-Coated Examination Gloves, XS	XS	Ivory	235 x 120 x 75	100	390 x 250 x 245	10
102646	Klinion Personal Protection Latex Polymer-Coated Examination Gloves, S	S	Ivory	235 x 120 x 75	100	390 x 250 x 245	10
102647	Klinion Personal Protection Latex Polymer-Coated Examination Gloves, M	M	Ivory	235 x 120 x 75	100	390 x 250 x 245	10
102648	Klinion Personal Protection Latex Polymer-Coated Examination Gloves, L	L	Ivory	235 x 120 x 75	100	390 x 250 x 245	10
102649	Klinion Personal Protection Latex Polymer-Coated Examination Gloves, XL	XL	Ivory	235 x 120 x 75	100	390 x 250 x 245	10

*Sizes are approximate

Product description

Disposable, powder-free, polymer-coated latex examination gloves

Intended use

The Klinion Protection examination gloves are intended to be worn on the hand for medical purposes and to help prevent contamination between patient and examiner.

User Info available at www.medeco.org

Classification

MDD 93/42/EEC Class I
PPE Regulation (EU) 2016/425 Category I
→ EU declaration of conformity at www.medeco.org

Features

- Economical and hygienic protection for hands
- Ambidextrous
- Textured surface
- Powder-free

MEDECO B.V.

Product Specifications

Product components

- Natural Rubber Latex
- Products listed on this specification are:
 - Phthalate-free
 - Absent of human or animal tissue
 - Absent of blood or blood derivative

Specifications

Dimensions (mm)

Width	XS	75 ± 5
	S	85 ± 5
	M	95 ± 5
	L	106 ± 5
	XL	116 ± 5
Thickness	Cuff	Min. 0.08
	Palm	Min. 0.10
	Finger	Min. 0.12
Length	Min. 240	

Physical Properties

Unaged	Tensile Strength	Min. 18 MPa
	Elongation at Break	Min. 650 %
	Force at Break	Min. 6.0 N
Aged	Tensile Strength	Min. 14 MPa
	Elongation at Break	Min. 500 %
	Force at Break	Min. 6.0 N
Freedom from Holes		AQL 1.5
Powder Residue		Max 2mg/glove
Protein		Max 50µg/g

MEDECO B.V.

Product Specifications

Quality Standards

- Manufactured in accordance with Quality Assurance ISO 13485 and ISO 9001

In compliance with:

 - EN 420:2003 + A1:2009 (fit for special purpose)
 - EN 455-1:2000, EN 455-2:2015, EN455-3:2015, EN455-4:2009

Sterilization

N/A – Product is provided non-sterile

Storage conditions

Store in a cool and dry place. Opened boxes should be kept away from fluorescent and sunlight. Gloves are packed in dispenser which is suitable for transport. Keep the gloves in the box when not in use.


Shelf life

3 years

Picture



Ref. product 102648

Approval	Name	Date	Signature
Product Manager	Levine van Buren	20-12-2017	
RA-Manager	Francis Dekker	20-12-2017	