

Paediatric Eco mask range

The lightest and most comfortable masks available



Intersurgical totally redesign the paediatric oxygen & aerosol mask

The lightest and most comfortable masks available

Every aspect of the design has been addressed, achieving unparalleled levels of comfort and fit. Not only comfortable for the patient but also comfortable for the environment.

The elimination of PVC from the masks has dramatically lowered their environmental impact.

Features and benefits

Incurved nose seal
improving the fit and eliminating the need for a separate nose clip. Designed to prevent oxygen or aerosolized drugs entering the patient's eyes

Soft face seal
for improved patient comfort

Offers a choice
of "under chin" or "on-chin" positions provides a better fit on a wider range of patient face shapes



Lightweight,
50% lighter than PVC equivalent

Greener 63% reduction in environmental impact

The shape and position of the exhalation ports have been designed to minimise the amount of rebreathed carbon dioxide

The material forming the body of the mask is clear, permitting visualisation of the patient's nose and lips. It is also rigid enough to maintain the mask's shape during transportation and use

Ordering information

Code	Description		Box qty
1190015	Paediatric Eco oxygen mask		50
1196015	Paediatric Eco oxygen mask with tube		50
1454015	Cirrus2 paediatric Eco mask kit & tube		36
1198015	Paediatric Eco aerosol mask		42
1138015	Paediatric Silente Eco mask kit		40

Completely PVC free options

1454010	Cirrus2 paediatric Eco mask kit with non-PVC tube		
1196010	Paediatric Eco oxygen mask kit with non-PVC tube		

IS10.27 • Issue 1 06.12