


Effective Pressure Area Care

HEEL PROTECTORS

High grade, antimicrobial TPU material with open design and ventilation holes to optimize microclimate. The Static Air® design ensures offloading of the heel and effective protection of the foot, ankle and lower limb.

Easy access to the heel, for dressing changes and visual inspection of the heel and foot plantar surface.

Silicone buttons to reduce infection risk and for easy cleaning.



4410 A Custom Heelprotector

Heel is continuously offloaded which complies to the latest EPUAP guidelines*.

Indications for use:

- to aid in the prevention of pressure injuries
- can be used for all levels of risk
- to aid in the treatment of all categories of pressure ulcers (I - IV)



LAD pump

EHOB® Heel Protectors have been clinically evaluated as part of effective protocol for the prevention and treatment of pressure injuries.

* EPUAP guidelines 2014 (<http://www.epuap.org/guidelines/quick-reference-guide-2014-edition-translations/>)



EHOB® has provided effective solutions for prevention and treatment of pressure injuries for more than 30 years.

We offer an extensive range of highly effective pressure relieving products to fit your specific needs.



EHOB® Heel Protectors are made of high grade, antimicrobial TPU material that can be cleaned and disinfected as per local policy guidelines.

Custom range available in 3 sizes for optimal fit for a wide range of patients. Standard range will fit most patients.



4450 A Standard Heel Protector with Velcro Straps and Cotton Lining



Optional with Light Air Delivery Pump (4450 B)

EHOB Static Air® Heel Protectors:

- Highly cost effective
- Easy to use
- Latex-free
- One year warranty
- Antimicrobial agent incorporated which is effective against a wide range of microorganisms including MRSA, C diff. and E coli..

Static Air® Heel Protector Range

Custom

- 4400 A Petite TPU unlined with Silicone Buttons
- 4410 A Standard TPU unlined with Silicone Buttons
- 4420 A XL TPU unlined with Silicone Buttons

Standard

- 4430 A One Size unlined TPU with Silicone Buttons
- 4440 A One Size lined TPU with Silicone Buttons
- 4450 A One Size lined TPU Velcro Straps

Add B (i.e. 4400 B) to order number to include LAD pump.
Also available: EHOB Truvue® Microfibre Heel Protector.
Please ask for the separate brochure.

Can be reused following local policy guidelines for cleaning and disinfection. All EHOB products in the Static Air range incorporate an antimicrobial agent that acts as an effective barrier against the following microorganisms: Methicillin-Resistant Staphylococcus Aureus (MRSA); Clostridium difficile; Escherichia coli; Klebsiella pneumonia; Pseudomonas aeruginosa; Salmonella choleraesuis; Candida albicans; Enterococcus faecium; Aspergillus niger and Bacillus subtilis.

