

Effective Pressure Area Care CUSHIONS

High grade, antimicrobial TPU material with optimal envelopment and immersion properties to ensure clinical efficacy and patient comfort.

Highly effective and complete pressure area care on any chair or wheelchair.

Indications for use:

- to aid in the prevention of pressure injuries
- can be used for all levels of risk
- to aid in the treatment of all categories of pressure ulcers (I - IV)

Ventilation Holes:


- aiding the even distribution of pressure across the patient's body surface area.
- provides effective ventilation and optimises micro-climate.
- minimises risk of skin maceration as body fluids are drawn away from the patient's skin.




4220 Pre-inflated Cushion (Opt. 4240 with Top-up Valve and Light Air Delivery Pump)

Pressure distribution results measured with EHOB Smart Seat System

Individual sitting on a chair with foam surface.



Same individual sitting on an EHOB Cushion on same chair.



EHOB® cushions have been clinically evaluated as part of effective protocol for the prevention and treatment of pressure injuries.



EHOB® has provided effective solutions for prevention and treatment of pressure injuries for more than 30 years.

We offer an extensive range of highly effective pressure relieving products to fit your specific needs.


EHOB® Cushions examples of use:



4260 TPU Multi Care Pad on Wheelchair



4230 TPU Standard Inflatable Cushion with LAD pump on Chair



Cushion Covers (optional):

4230 COV for 4210/4230
4240 COV for 4220/4240
4250 COV for 4250

EHOB® Cushions:

- Highly cost effective
- Easy to use
- Latex free
- Low profile (2 to 5 cm)
- One year warranty
- Up to 320 kg weight limit
- Antimicrobial agent incorporated which is effective against a wide range of microorganisms including MRSA, C. Diff and E.Coli



EHOB® Cushion Range

Cushions

4210A	43x43x5 cm	136 Kg	TPU Standard Cushion pre-inflated
4220A	51x51x5 cm	136 Kg	TPU Large Cushion pre-inflated
4230A	43x43x4 cm	136 Kg	TPU Standard Cushion with valve
4230B	43x43x4 cm	136 Kg	TPU Standard Cushion with valve & pump
4240A	51x51x5 cm	136 Kg	TPU Large Cushion with valve
4240B	51x51x5 cm	136 Kg	TPU Large Cushion with valve & pump
4250	71x56x6 cm	320 Kg	TPU Bariatric Cushion with valve & pump
4260	91x38 cm	136 Kg	TPU Multi Care Pad with valve & pump

Straps

4210Z	Cushion fixture strap for 4210
4220Z	Cushion fixture strap for 4220
4220ZA	Cushion fixture straps for 4210/4220 6-pack

Can be reused following local policy guidelines for cleaning and disinfection. All EHOB products in the Static Air range incorporate an antimicrobial agent that acts as an effective barrier against the following microorganisms: Methicillin-Resistant Staphylococcus Aureus (MRSA); Clostridium difficile; Escherichia coli; Klebsiella pneumonia; Pseudomonas aeruginosa; Salmonella choleraesuis; Candida albicans; Enterococcus faecium; Aspergillus niger and Bacillus subtilis.

EHOB®
prevention • treatment • solutions



EHOB Europe - Lichtenauerlaan 102-120
3062 ME Rotterdam
T: +31 (0)10 7892404
E: paul.mersel@ehob.com

ehob-europe.com