


ASSORTIMENTSOVERZICHT TERUMO

INJECTIENAALDEN

	Art.nr	Afmeting (mm)	Afmeting (gauge)	ZI-nummer	Verpakking
	9080150	0,4 x 16 mm	27G x 5/8	15700240	100 st.
	9080125	0,4 x 20 mm	27G x 3/4	15691616	100 st.
	9080124	0,45 x 12 mm	26G x 1/2	15700224	100 st.
	9080122	0,5 x 16 mm	25G x 5/8	15691357	100 st.
	9080123	0,5 x 25 mm	25G x 1	15700216	100 st.
	9080119	0,6 x 25 mm	23G x 1	15691365	100 st.
	9080116	0,7 x 30 mm	22G x 1 1/4	15691403	100 st.
	9080112	0,8 x 16 mm	21G x 5/8	15691446	100 st.
	9080114	0,8 x 38 mm	21G x 1 1/2	15691462	100 st.
	9080115	0,8 x 50 mm	21G x 2	15691470	100 st.
	9080106	0,9 x 25 mm	20G x 1	15700143	100 st.
	9080108	0,9 x 40 mm	20G x 1 1/2	15691489	100 st.
	9080109	0,9 x 50 mm	20G x 2	15691500	100 st.
	9080104	1,1 x 40 mm	19G x 1 1/2	15700127	100 st.
	9080100	1,2 x 40 mm	18G x 1 1/2	15691527	100 st.


INJECTIENAALDEN

INTRACUTAAN (in de huid)




Maat (mm)	Kleur	ZI-nr.	Art.nr
0,4 x 16	grijs	15700240	9080150
0,45 x 12	bruin	15700224	9080124
0,5 x 16	oranje	15691357	9080122
0,8 x 16	groen	15691446	9080112

SUBCUTAAN (in het bindweefsel)




Maat (mm)	Kleur	ZI-nr.	Art.nr
0,4 x 16	grijs	15700240	9080150
0,4 x 20	grijs	15691616	9080125
0,45 x 12	bruin	15700224	9080124
0,5 x 16	oranje	15691357	9080122
0,5 x 25	oranje	15700216	9080123
0,6 x 25	blauw	15691365	9080119
0,8 x 16	groen	15691446	9080112
0,9 x 25	geel	15700143	9080106

INTRAMUSCULAIR (in het spierweefsel)



Maat (mm)	Kleur	ZI-nr.	Art.nr
0,4 x 20	grijs	15691616	9080125
0,5 x 25	oranje	15700216	9080123
0,6 x 25	blauw	15691365	9080119
0,7 x 30	zwart	15691403	9080116
0,8 x 38	groen	15691462	9080114
0,8 x 50	groen	15691470	9080115
0,9 x 25	geel	15700143	9080106
0,9 x 40	geel	15691489	9080108
0,9 x 50	geel	15691500	9080109
1,1 x 40	crème	15700127	9080104
1,2 x 40	rose	15691527	9080100

INTRAVENEUS (in de aderen)



Maat (mm)	Kleur	ZI-nr.	Art.nr
0,6 x 25	blauw	15691365	9080119
0,7 x 30	zwart	15691403	9080116
0,8 x 38	groen	15691462	9080114
0,9 x 25	geel	15700143	9080106
0,9 x 40	geel	15691489	9080108
1,1 x 40	crème	15700127	9080104
1,2 x 40	rose	15691527	9080100

